

EVERY CONNECTION COUNTS



For 45 years, Baxter has offered a comprehensive line of syringes designed to **prevent misconnection** with other devices. In an effort to prevent wrong route delivery (misconnections) of fluids and gases, the International Organization for Standardization (ISO), in collaboration with the FDA, AAMI, clinicians, has developed a new standard to reduce the likelihood of tubing misconnections, ISO 80369. Baxter, in support of this global initiative to **enhance patient safety**, is introducing a complete line of **ENFit** enteral syringes and accessories that complies with ISO 80369-3 specific for enteral devices.

To comply with ISO 80369-3 and help ensure patient safety, GEDSA is recommending to manufacturers and healthcare providers the adoption of new **ENFit** connectors in North America.

To facilitate conversion to the **ENFit** system, Baxter is offering a **transition connector** that allows clinicians to use the new administration and patient tube sets with the existing male tip syringes. Baxter's transition connector is available in non-sterile and sterile configurations (cap included but not shown).

Contact Baxter Medical Information at 1-800-422-2751 or your Baxter Sales Representative for more information.

ORDERING INFORMATION



REF

ENFit Enteral-to-Oral Connector ¹	30	H93842907EN
Sterile ENFit Enteral-to-Oral Connector ²	30	H93842905EN

1 Available in North America. Coming soon to Australia and New Zealand.

2 Available for Western Europe.

For more information about ISO 80369-3 and the Stay Connected Initiative, visit www.stayconnected.org.

www.baxter.com

Baxter Healthcare Corporation
One Baxter Parkway
Deerfield, IL 60015

Baxter is a trademarks of Baxter International Inc., or its subsidiaries.
ENFit is a registered trademark of Global Enteral Device Supplier Association, Inc.

USMP/MG100/16-0002 12/16

MAKE ORAL CONNECTIONS

We've created an Oral Liquid Medication Delivery portfolio with a full range of syringes, sizes, packaging configurations and accessories so that you have access to the appropriate drug delivery system that works best for your pharmacy needs. Our goal is to help you streamline your pharmacy practices. That is why we developed a complete system for preparing and administering oral medications with safety and efficiency in mind.



Baxter